
ULT JUMBO Filtertrolley 2.0 MED

ULT JUMBO®

Technical Documentation for Extraction and Filter Unit **Type Series Jumbo MED**

During medical applications like laser surgery, electrocautery and others hazardous surgical plume, viruses and bacterias in the form of dust, floating particles, gases and fumes are generated. This smoke can be a respiratory irritant to the surgical team and to the patients, but more it can be a source of airborne infection. To minimize the risk of infectious particulates and also according to regulations of air purification at the workplace it is recommended that a smoke evacuator is used.

The smoke evacuator ULT JUMBO Filtertrolley is especially designed for these tasks. The optional recirculation ventilation helps to reduce energy costs. The simple and application-oriented, partially even individual, solutions result from longtime experience in theory and practice and a multitude of different application cases.

Advantages of the JUMBO Filtertrolley

- Ultra quiet
- compact design
- high performance fans with low energy consumption
- user-friendly and cost-effective
- easy operation, maintenance and filter replacement

Fields of Application

Medical laser surgery and electrocautery , Laboratories
Sterilization rooms, Patient's awakening rooms
Anaesthesia rooms, Dentist's surgery

Possible Connections

Extraction hoods and tubes



Filter Technology

- prefilter carton F 5 – disposal is free from contamination
- HEPA filter H 13
- Activated carbon filter

Technical Documentation for Extraction and Filter Unit Type Series Jumbo MED

1. inlet stub
2. filter cartridge
3. air outlet
4. motor / fan
5. flow rate controller
6. power switch
7. loaded filter indicator

Functional principle:

A special fan (4) with a large capacity of spare pressure creates an air flow adapted to the respective application at the clean air section of the filter cartridge (2). An optional flow rate controller is available. The hazardous air is reliably extracted (1) and cleaned. Because of the high purity of the filtered air it can be supplied back to the working area (3), thus, reducing energy costs.

The different filters are combined in one filter cartridge and can be replaced singularly. The filter is recyclable.


Technical Documentation for Extraction and Filter Unit Type Series Jumbo MED

Specification / Fan Graph

Type Series		Jumbo
Parameters	Unit	
max. flow rate	m ³ /h	170
max. vacuum	Pa	2.800
nominal capacity	m ³ /h at Pa	80 / 1.400
nominal motor power	kW	0,15
nominal voltage	V	230
nominal current	A	0,95
frequency	Hz	50
type of protection	IP	54
noise level at 50%	dB(A)	45
noise level at 100%	dB(A)	49
weight	kg	18
flow rate controller		optional
loaded filter indicator		optional
inlet stub	mm	1 x Ø50 or 1 x Ø75 or 2 x Ø50
air outlet		2 x adjustable nozzles

Type series Jumbo Filtertrolley 2.0


Technical Documentation for Extraction and Filter Unit Type Series Jumbo MED

Dimensions Filtertrolley 2.0 MED


alternatively: 1x Ø 50mm
1x Ø 75mm
2x Ø 50mm


Technical Documentation for Extraction and Filter Unit Type Series Jumbo MED

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Brevo B.V.